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## 1. PACKAGE LEAFLET

### URGO Mouth ulcers & wounds

Spray indicated to help treat mouth ulcers, gingivitis, burns and wounds related to dental procedures scaling and braces.

Ideal for both adults and children (can be used from 36 months), this product is sugar and alcohol-free.

It does not sting.

### How does URGO Mouth ulcers & wounds work?

Containing low-molecular weight hyaluronic acid (0.5%), URGO Mouth ulcers & wounds is a medical device that protects hydrates and promotes repair of the mucous membranes of the mouth in the event of mouth ulcers or wounds. As a result it relieves pain and promotes healing.

### Instructions for use:

When using for the first time, prime the pump mechanism by pressing once or twice. Point the nozzle towards the mouth wound and spray 2 to 3 times in order to form a thin, even layer. Repeat at least 3 times per day to ensure continued protection of the wound. For optimal efficacy, use the spray after meals and after tooth-brushing.

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**Kontraindikace:**

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**Precautions for use:**

Rinse the nozzle with warm water, dry it with a clean paper and fold back the nozzle after each use.

Stop treatment in the event of irritation or hypersensitivity to the product and consult a doctor if symptoms persist.

Keep out of the reach of children.

Use in children must be supervised by an adult.

Check that the packaging and bottle are intact before use.

Store at a temperature of between 5°C and 40°C.

**Contraindications:**

Children aged under 3 years.